<p>EN</p> <p>Verified body responsible for certification and ongoing conformity:</p> <p>SATRA Technology Europe Ltd Bracetown Business Park Clonee, Dublin D15 YN2P, Ireland (CE 2777)</p>	<p> Dongguan Grinvald Technology Co. Ltd. 401, Building #3, No 4 of Guangming New Village 2 Road Dongcheng Dongguan City Guangdong China 523000</p>	<p> Global Tooling Service s.r.o. Na Follmance 2155/15, 120 00 Praha 2 - Vinohrady, Czech Republic</p>
--	---	--

Product reference: MD0120	Sizes available: S(6), M(7), L(8), XL(9)
----------------------------------	---

INSTRUCTIONS FOR USE

Declaration:

The manufacturer declares that the properties of the given medical device and personal protective equipment meet all the requirements pursuant to Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and Regulation (EU) 2016/425 on personal protective equipment.

Intended use:

Nitrile medical examination gloves are non-sterile and are used to perform medical examinations, diagnostic or therapeutic procedures, and to work with contagious medical materials. They protect the patient and user from contamination. The gloves are intended for single use. The gloves can also be used for contact with food.

Classification of gloves pursuant to legislation:

Risk class of the medical device: I
Personal protective equipment risk category: III

Harmonised standards used: EN ISO 13485:2016, EN ISO 9001:2015, EN ISO 14971:2019, IEC 62366-1:2015(en), EN ISO 15223-1:2016, EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, EN 16523-1:2015+A1:2018, EN ISO 374-1:2016+A1:2018, EN 420:2003+A1:2009, EN 374-2:2014, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 1186:2020

Results and Restrictions on Use - This product has been tested with the following results:

EN ISO 374-1:2016+A1:2018 /Typ C	Level	EN ISO 374-4:2019 Degradation %	EN ISO 374-1/Type C
40% sodium hydroxide (K)	6	- 18.0	 
EN ISO 374-5:2016 Protection against bacteria and fungi Virus protection	Pass Pass		EN ISO 374-5:2016   VIRUS

EN ISO 374-1:2016+A1:2018
Permeation levels are based on penetration times as follows:

Permeation level	1	2	3	4	5	6
Measured penetration time (min.)	>10	>30	>60	>120	>240	>480

EN ISO 374-4:2019 The degradation results indicate a change in the gloves' puncture resistance after exposure to a strong chemical:

EN ISO 374-5:2016 Penetration resistance was assessed under laboratory conditions and applies to the test samples only.

Warning:

This information does not reflect the actual period of protection in the workplace and the differences between mixtures and pure chemicals.

Chemical resistance was tested under laboratory conditions and only on samples taken from the palm (except when the gloves are equal to or exceed 400 mm - when the cuffs were also tested) and only applies to the test chemical. If a chemical is used in a mixture, the result may be different.

We recommend checking that the gloves are suitable for the intended use. Workplace conditions may differ from those at the type test depending on temperature, abrasion and degradation.

When in use, the gloves may be less resistant to hazardous chemicals due to changes in physical properties. Movements, snagging, friction, degradation due to chemical contact can significantly shorten the actual period of use/lifespan. When using corrosive chemicals, degradation can be the most important aspect for choosing chemically resistant gloves.

Additional information:

Inspect before use, do not use damaged gloves. The gloves are disposable and ambidextrous. Protect from ozone. Use by the date shown on the product packaging. Store in a dry place and keep away from sunlight. Storage temperature 10 °C - 30 °C.

The gloves do not contain any substances known to cause allergies. The gloves do not provide any mechanical protection. Disposable, do not throw in waste bin. Check for damage before use, do not use damaged gloves.

The declaration of conformity can be found on the website: <https://medcaregloves.com/>

Putting on:

- Before putting on, remove all jewellery from your hands and wrists and wash your hands.
- Place the gloves on a prepared work surface.
- Place the glove on your dominant hand by holding the glove with your other hand, remembering to touch only the inside of the glove and pulling it onto the dominant hand until you reach the level of your fingers.
- Use the dominant gloved hand to put the other glove on the non-dominant hand.
- Once both gloves are on, you can touch the outside of the gloves to adjust their fit.

Taking off:

- Using the dominant hand, start by holding the outside of the glove on the non-dominant hand by the side of the palm near the cuff.
- Remove the glove from the non-dominant hand, grasp it in the gloved hand and squeeze it.
- Insert two fingers under the glove cuff on the other hand and carefully pull it away from the hand without touching the wrist. Pull this glove so that it turns inside out while wrapping the first glove inside it.
- Dispose of the gloves in a designated waste container.

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